General name	Electric therapy equipment	El
Brand Name	RELTEC MD21	De
Model	MD21	pr
External dimensions	300(W)×110(H)×280(D)(mm)	D
Mains voltage	100-240 V ~	0
Rated frequency	50-60 Hz	
Power consumption	50VA	
Weight	Approx. 2.2Kg (main body)	(
Treatment time Manual setting (1 to 99 minutes)		
Output voltage	DC-5500±500 V	(
Essential Performance	The essential performance of this product is "output voltage DC-5500V \pm 500V". If the essential performance is lost or deteriorated due to electromagnetic interference, the display may flash or become malfunction. In this case, turn off the power of the device and stop the treatment.	
Output pole	Available options: monopole mode, dipole mode	
Safety features	The device automatically stops when a conductive plate plug is disconnected. The device automatically stops when a conductive plate plug is not properly inserted. The device automatically stops when the main body is installed improperly.	

	classification		Class I	
	Degree of electric shock protection		BF-type applied part	
	Degree of protection		IPX0	
	Operation mode		Consequtive operation	
	Usage /Storage environment	operating temperature	10~39℃	
		operating humidity	30∼85%RH (no condensation))	
		transport and storage temperature	0~50℃	
	onment	transport and storage humidity	10~90%RH (no condensation)	

tric shock protection

70	Type BD	Type YD
Accessories	Conductive plate AC (2.5m) 2 pcs Conductive plate BD (2.5m) 2 pcs Cover (S) 4 pcs Cover (XS) 4 pcs Velcro XS (30-cm length) 2 pcs Velcro M (80-cm length) 1 pc Velcro L (120-cm length) 1 pc AC adapter (1.25m) 1 pc Power cord (2 m) 1 pc Instruction manual 1 copy Insulating mat (90 cm×190 cm) 1 pc Spiral tube 10pc	Conductive plate AC (2.5m) 2 pcs Conductive plate YD (2.5m) 2 pcs Cover (S) 4 pcs Cover (XS) 8 pcs Velcro XS (30-cm length) 1 pc Velcro M (80-cm length) 1 pc Velcro L (120-cm length) 1 pc AC adapter (1.25m) 1 pc Power cord (2 m) 1 pc Instruction manual 1 copy Insulating mat (90 cm×190 cm) 1 pc ENERGY CHECK 1 pc Spiral tube 10.5m) 10pc

User Profile

Operator	Doctor, Nurse , Physiotherapist Occupational therapist	Patient	
Gender	Men and Women	Gender	≧3 years old
Age	≧ 20 years old	Weight	≧10 kg
Professional Competence	Medical License, Nurse's License, Physiotherapy certificate, Occupational therapist certificate	Part of the body applied to	Limb, Shoulder, Back of head, Abdomen
Use Frequency	30 minutes/time several times/ day	Location of Use	Professional healthcare facility environment
Linguistic	Japanese	Location of Ose	
Knowledge	At least 11 years old – 5 years intensive reading experience(school)		